



Should I buy the NI card?

No. You can purchase the NI card on your own. In this case all the guarantee, delivery and installation duties are on you.

Can I use more than one software on a PC?

Yes. All the members of the SPEL software family are developed on the same type of NI card.

Can I use laptop?

Yes. The PCMCIA version of the NI card is available.

Do I have upgrade?

Maintenance upgrades are available. Further developed versions are treated as new products. These are available for a moderate upgrade cost.

Can I have software training?

Yes. For detailed information please contact Experimetria.

Do I have guarantee?

Yes. 2 years.

What is in the box?

The box contains:

- SPEL software
- SPEL software user manual
- NI Daq A/D card
- NI Daq driver software
- AIF-0# connector box
- cables

Is it an open system?

Can I use my biological amplifiers?

Yes. The SPEL software systems can receive any standard analogue 0 +5v DC signal.

Can I measure other parameters?

Yes. All the SPEL software systems contains the User Define option, which is capable to receive any custom define voltage data.

Who gives the computer?

The Experimetria do not supply the PC with the software system instantly. In case of request Experimetria kindly delivers the PC or laptop.

Can I use the software on an other PC or laptop?

No. If you would like to analyse the previously recorded data Experimetria kindly provide you the demo version of your software where the full analysis option is presented.

What is the Windows operation system?

The SPEL software systems can work under Win 9X, NT, ME and XP. The Win 2000 is not applicable.